



**The Streamlight EMS ProTac provides the optimal non-glare light for pupil and wound examination.**

|  |  |
|--|--|
| <p>Powers on at low beam: 3.6 lumens/runs 36 hours; ideal for pupil exams. Faint peripheral illumination surrounding the focused beam slightly illuminates the opposite pupil without directly causing it to contract, facilitating observation of consensual contraction in dim light</p> | <p>One-handed operation of tailcap switch; easy to use with gloved hands</p> |
| <p>Two taps of the tail switch accesses medium mode, which delivers 10 lumens/runs for 19 hours. Ideal for non-glare wound examination</p>   | <p>IPX7 rated design; waterproof to 1 meter for 30 minutes</p>               |
| <p>Three taps of the tail switch accesses high mode, which delivers 50 lumens/runs 1.5 hours (alkaline); runs 3.4 hours using AA lithium battery</p>   | <p>Removable pocket clip; includes nylon holster</p>                         |
| <p>LED Solid State power regulation provides maximum light output throughout battery life</p>  | <p>Designed by a medical services professional</p>                           |
|  | <p>3.97" (10.08 cm); 2.3 oz (65 grams)</p>                                   |
|  | <p><b>CE</b> Meets applicable European Community Directives</p>              |
|  | <p>Limited lifetime warranty</p>   |
|  | <p>#88034 - EMS ProTac with 1 "AA" alkaline battery and nylon holster</p>    |